

A18133

Leader in Biomolecular Solutions for Life Science



ANGPTL8 Rabbit pAb

Catalog No.: A18133

Basic Information

Observed MW

24kDa

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to enable hormone activity. Involved in regulation of lipid metabolic process and triglyceride homeostasis. Acts upstream of or within positive regulation of protein processing and regulation of lipoprotein metabolic process. Located in extracellular region.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

55908

Swiss Prot

Q6UXH0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 22-198 of human ANGPTL8 (NP_061157.3).

Synonyms

RIFL; TD26; PRO1185; PVPA599; C19orf80; ANGPTL8

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates, using ANGPTL8 Rabbit pAb (A18133) at 1:3000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 180s.