

A18193

Leader in Biomolecular Solutions for Life Science



CFAP69 Rabbit pAb

Catalog No.: A18193

Basic Information

Observed MW

Refer to figures

Calculated MW

106kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Acts upstream of or within sperm axoneme assembly. Located in cytoplasm and sperm midpiece. Implicated in spermatogenic failure 24.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

79846

Swiss Prot

A5D8W1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

FAP69; SPGF24; C7orf63; CFAP69

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.