

A18193

Leader in Biomolecular Solutions for Life Science



CFAP69 Rabbit pAb

Catalog No.: A18193

Basic Information

Observed MW

Refer to figures

Calculated MW

106kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Acts upstream of or within sperm axoneme assembly. Located in cytoplasm and sperm midpiece. Implicated in spermatogenic failure 24.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79846

Swiss Prot

A5D8W1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human CFAP69 (NP_001034795.2).

Synonyms

FAP69; SPGF24; C7orf63; CFAP69

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.