

A18230

Leader in Biomolecular Solutions for Life Science



CLTCL1 Rabbit pAb

Catalog No.: A18230

Basic Information

Observed MW

Refer to figures

Calculated MW

187kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene is a member of the clathrin heavy chain family and encodes a major protein of the polyhedral coat of coated pits and vesicles. Chromosomal aberrations involving this gene are associated with meningioma, DiGeorge syndrome, and velo-cardio-facial syndrome. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

8218

Swiss Prot

P53675

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 1600-1663 of human CLTCL1 (NP_009029.3).

Synonyms

CLTD; CHC22; CLH22; CLTCL; CLTCL1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.