

A18413

Leader in Biomolecular Solutions for Life Science



MAGEC1 Rabbit pAb

Catalog No.: A18413

Basic Information

Observed MW

Refer to figures

Calculated MW

124kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

This gene is a member of the melanoma antigen gene (MAGE) family. The proteins of this family are tumor-specific antigens that can be recognized by autologous cytolytic T lymphocytes. This protein contains a large number of unique short repetitive sequences in front of the MAGE-homologous sequence, and therefore is about 800 aa longer than the other MAGE proteins.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9947

Swiss Prot

O60732

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 600-700 of human MAGEC1 (NP_005453.2).

Synonyms

CT7; CT7.1; MAGEC1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.