

A18422

Leader in Biomolecular Solutions for Life Science



ANKRD26 Rabbit pAb

Catalog No.: A18422

Basic Information

Observed MW

Refer to figures

Calculated MW

196kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

This gene encodes a protein containing N-terminal ankyrin repeats which function in protein-protein interactions. Mutations in this gene are associated with autosomal dominant thrombocytopenia-2. Pseudogenes of this gene are found on chromosome 7, 10, 13 and 16. Multiple transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
22852

Swiss Prot
Q9UPS8

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

THC2; bA145E8.1; ANKRD26

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.