

A18435

Leader in Biomolecular Solutions for Life Science



# ANKS1A Rabbit pAb

Catalog No.: A18435

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

123kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human

## Background

---

Predicted to enable ephrin receptor binding activity. Predicted to be involved in ephrin receptor signaling pathway; neuron remodeling; and substrate-dependent cell migration. Predicted to act upstream of or within negative regulation of ubiquitin-dependent protein catabolic process and regulation of ephrin receptor signaling pathway. Located in cytosol and nucleoplasm.

## Recommended Dilutions

---

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

---

### Gene ID

23294

### Swiss Prot

Q92625

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

ANKS1; ANKS1A

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.