

A18477

Leader in Biomolecular Solutions for Life Science



NAT14 Rabbit pAb

Catalog No.: A18477

Basic Information

Observed MW

Refer to figures

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable N-acetyltransferase activity. Predicted to be integral component of membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

57106

Swiss Prot

Q8WUY8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 77-206 of human NAT14 (NP_065111.1).

Synonyms

KLP1; NAT14

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.