

A18515

Leader in Biomolecular Solutions for Life Science



NAA11 Rabbit pAb

Catalog No.: A18515

Basic Information

Observed MW

Refer to figures

Calculated MW

26kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Enables peptide alpha-N-acetyltransferase activity. Involved in N-terminal protein amino acid acetylation. Located in several cellular components, including Golgi apparatus; centrosome; and nucleoplasm. Part of NatA complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

84779

Swiss Prot

Q9BSU3

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 170-229 of human NAA11 (NP_116082.1).

Synonyms

ARD2; ARD1B; hARD2; NAA11

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.