

A18523

Leader in Biomolecular Solutions for Life Science



ZNF276 Rabbit pAb

Catalog No.: A18523

Basic Information

Observed MW

Refer to figures

Calculated MW

67kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse,Rat

Background

Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Located in kinetochore.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

92822

Swiss Prot

Q8N554

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 180-360 of human ZNF276 (NP_689500.2).

Synonyms

ZADT; CENPZ; CENP-Z; ZFP276; ZNF477; ZNF276

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.