

A18523

Leader in Biomolecular Solutions for Life Science



ZNF276 Rabbit pAb

Catalog No.: A18523

Basic Information

Observed MW

Refer to figures

Calculated MW

67kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse,Rat

Background

Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Located in kinetochore.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
92822

Swiss Prot
Q8N554

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ZADT; CENPZ; CENP-Z; ZFP276; ZNF477; ZNF276

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.