

A18535

Leader in Biomolecular Solutions for Life Science



# ACOT4 Rabbit pAb

Catalog No.: A18535

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

46kDa

### Category

Polyclonal Antibody

### Applications

IHC-P,ELISA

### Cross-Reactivity

Human,Mouse,Rat

## Background

Enables acyl-CoA hydrolase activity and succinyl-CoA hydrolase activity. Involved in carboxylic acid metabolic process; saturated monocarboxylic acid metabolic process; and succinyl-CoA metabolic process. Located in peroxisome.

## Recommended Dilutions

IHC-P 1:50 - 1:200

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

122970

### Swiss Prot

Q8N9L9

### Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

### Synonyms

PTE1B; PTE2B; PTE-Ib; ACOT4

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.