

A18548

Leader in Biomolecular Solutions for Life Science



CD300LG Rabbit pAb

Catalog No.: A18548

Basic Information

Observed MW

Refer to figures

Calculated MW

36kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

Members of the CD300 (see MIM 606786)-like (CD300L) family, such as CD300LG, are widely expressed on hematopoietic cells. All CD300L proteins are type I cell surface glycoproteins that contain a single immunoglobulin (Ig) V-like domain (Takatsu et al., 2006 [PubMed 16876123]).

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

146894

Swiss Prot

Q6UXG3

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

CLM9; CLM-9; TREM4; TREM-4; NEPMUCIN; CD300LG

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.