

A18556

Leader in Biomolecular Solutions for Life Science



ZCCHC12 Rabbit pAb

Catalog No.: A18556

Basic Information

Observed MW

Refer to figures

Calculated MW

45kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a downstream effector of bone morphogenetic protein (BMP) signalling. This protein contains a zinc finger domain and functions as a transcriptional coactivator. Variation in this gene may be associated with X-linked cognitive disability. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

170261

Swiss Prot

Q6PEW1

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

SIZN; SIZN1; PNMA7A; ZCCHC12

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.