

A18568

Leader in Biomolecular Solutions for Life Science



TAZ1 Rabbit pAb

Catalog No.: A18568

Basic Information

Observed MW

191kDa

Calculated MW

191kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Predicted to be involved in spermatid development and spermatocyte division. Predicted to act upstream of or within apoptotic process; ncRNA transcription; and positive regulation of meiotic cell cycle phase transition. Predicted to be located in cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

375337

Swiss Prot

Q8N9V7

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 13-276 of mouse TAZ1 (NP_001186665.1).

Synonyms

C3orf77; TAZ1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from Mouse brain, using TAZ1 Rabbit pAb (A18568) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.