

A19211

Leader in Biomolecular Solutions for Life Science



GPA5 Rabbit pAb

Catalog No.: A19211

Basic Information

Observed MW

120kDa

Calculated MW

80kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Oryza sativa

Background

GPA5 is a membrane-associated protein capable of forming homodimers and that it is specifically localized to DVs in developing endosperm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

4341636

Swiss Prot

Q67WV5

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Rab5a Effector; GPA5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of various tissues from the japonica rice (*Oryza sativa* L.) variety Zhonghua 11, using GPA5 Rabbit pAb (A19211) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Enhanced Kit (RM00021).

Exposure time: 30s.