

A19889

Leader in Biomolecular Solutions for Life Science



[KO Validated] TXLNA Rabbit pAb

Catalog No.: A19889 **KO Validated**

Basic Information

Observed MW

80kDa

Calculated MW

62kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to enable syntaxin binding activity. Predicted to be involved in exocytosis. Predicted to act upstream of or within B cell activation. Located in cytoplasm.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

200081

Swiss Prot

P40222

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

IL14; TXLN; NA

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of lysates from wild type (WT) and TXLNA knockout (KO) HeLa cells, using [KO Validated] TXLNA Rabbit pAb (A19889) at 1:1000 dilution. Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.