

A20005

Leader in Biomolecular Solutions for Life Science



[KO Validated] USP4 Rabbit pAb

Catalog No.: A20005 **KO Validated**

Basic Information

Observed MW

108kDa

Calculated MW

109kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Mouse

Background

The protein encoded by this gene is a protease that deubiquitinates target proteins such as ADORA2A and TRIM21. The encoded protein shuttles between the nucleus and cytoplasm and is involved in maintaining operational fidelity in the endoplasmic reticulum. Three transcript variants encoding different isoforms have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

7375

Swiss Prot

Q13107

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 560-700 of human USP4 (NP_003354.2).

Synonyms

UNP; Unph; P4

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from wild type (WT) and USP4 knockout (KO) HeLa cells, using [KO Validated] USP4 Rabbit pAb (A20005) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 1s.