

A20291

Leader in Biomolecular Solutions for Life Science



ZW10 Rabbit pAb

Catalog No.: A20291

Basic Information

Observed MW

Refer to figures

Calculated MW

89kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a protein that is one of many involved in mechanisms to ensure proper chromosome segregation during cell division. This protein is an essential component of the mitotic checkpoint, which prevents cells from prematurely exiting mitosis.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9183

Swiss Prot

O43264

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HZW10; KNTC1AP; ZW10

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.