

A20339

Leader in Biomolecular Solutions for Life Science



ASH1L Rabbit pAb

Catalog No.: A20339

Basic Information

Observed MW

Refer to figures

Calculated MW

333kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of the trithorax group of transcriptional activators. The protein contains four AT hooks, a SET domain, a PHD-finger motif, and a bromodomain. It is localized to many small speckles in the nucleus, and also to cell-cell tight junctions.

Recommended Dilutions

WB 1:1000 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
55870

Swiss Prot
Q9NR48

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

ASH1; KMT2H; MRD52; ASH1L1; ASH1L

Contact

 www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.