

A20361

Leader in Biomolecular Solutions for Life Science



RASSF7 Rabbit pAb

Catalog No.: A20361

Basic Information

Observed MW

40kDa

Calculated MW

40kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human,Rat

Background

Predicted to be involved in apoptotic process; regulation of microtubule cytoskeleton organization; and signal transduction. Located in centriolar satellite.

Recommended Dilutions

WB 1:500 - 1:1000

Elisa 1:500 - 1:2000

Immunogen Information

Gene ID

8045

Swiss Prot

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HRC1; FAP88; HRAS1; CFAP88; C11orf13; RASSF7

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.