

A20430

Leader in Biomolecular Solutions for Life Science



PRSS35 Rabbit pAb

Catalog No.: A20430

Basic Information

Observed MW

Refer to figures

Calculated MW

47kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to be located in extracellular region.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

167681

Swiss Prot

Q8N3Z0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

C6orf158; dJ223E3.1; PRSS35

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.