

A20648

Leader in Biomolecular Solutions for Life Science



STK17B Rabbit pAb

Catalog No.: A20648

Basic Information

Observed MW

Refer to figures

Calculated MW

42kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables ATP binding activity and protein serine/threonine kinase activity. Involved in intracellular signal transduction; positive regulation of fibroblast apoptotic process; and protein phosphorylation. Located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

9262

Swiss Prot

O94768

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

DRAK2; STK17B

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.