

A20780

Leader in Biomolecular Solutions for Life Science



ADGRG3 Rabbit pAb

Catalog No.: A20780

Basic Information

Observed MW

Refer to figures

Calculated MW

61kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Predicted to enable G protein-coupled receptor activity. Involved in G protein-coupled receptor signaling pathway and regulation of cell migration. Located in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

222487

Swiss Prot

Q86Y34

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 21-270 of human ADGRG3 (NP_740746.4).

Synonyms

PB99; GPR97; PGR26; ADGRG3

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.