

A2119

Leader in Biomolecular Solutions for Life Science



# UAP1 Rabbit pAb

Catalog No.: A2119

## Basic Information

### Observed MW

### Calculated MW

59kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Rat

## Background

Enables identical protein binding activity. Predicted to be involved in UDP-N-acetylglucosamine biosynthetic process. Located in cytosol; nucleoplasm; and plasma membrane.

## Recommended Dilutions

**WB** 1:500 - 1:2000

**ELISA** Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

6675

### Swiss Prot

Q16222

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human UAP1 (NP\_003106.3).

### Synonyms

AGX; AGX1; AGX2; SPAG2; UAP1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.