

A21464

Leader in Biomolecular Solutions for Life Science



IKK β Rabbit pAb

Catalog No.: A21464

1 Publications

Basic Information

Observed MW

87kDa

Calculated MW

87kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Mouse

Background

The protein encoded by this gene phosphorylates the inhibitor in the inhibitor/NF-kappa-B complex, causing dissociation of the inhibitor and activation of NF-kappa-B. The encoded protein itself is found in a complex of proteins. Several transcript variants, some protein-coding and some not, have been found for this gene.

Recommended Dilutions

WB 1:500 - 1:1000

ELISA Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
3551

Swiss Prot
O14920

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

IKK2; IKKB; IMD15; IMD15A; IMD15B; NFKB1KB; IKK-beta; IKK β

Contact



www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.