

A21465

Leader in Biomolecular Solutions for Life Science



[KO Validated] INPPL1 Rabbit pAb

Catalog No.: A21465

KO Validated

Basic Information

Observed MW

Refer to figures

Calculated MW

139kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

The protein encoded by this gene is an SH2-containing 5'-inositol phosphatase that is involved in the regulation of insulin function. The encoded protein also plays a role in the regulation of epidermal growth factor receptor turnover and actin remodelling. Additionally, this gene supports metastatic growth in breast cancer and is a valuable biomarker for breast cancer.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3636

Swiss Prot

O15357

Immunogen

Recombinant protein of human INPPL1.

Synonyms

OPSM; SHIP2; L1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.