

A22102

Leader in Biomolecular Solutions for Life Science



MEX3A Rabbit pAb

Catalog No.: A22102 **1 Publications**

Basic Information

Observed MW

54kDa

Calculated MW

57kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Human

Background

Enables RNA binding activity. Located in P-body and cytosol.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

92312

Swiss Prot

P08684

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

RKHD4; MEX-3A; RNF162; MEX3A

Contact

 www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.