

A22795

Leader in Biomolecular Solutions for Life Science



# PE Mouse anti-Human CD3 mAb

Catalog No.: A22795

1 Publications

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

23kDa

### Category

Monoclonal Antibody

### Applications

FC

### Cross-Reactivity

Human

### CloneNo number

AMC0533

### Conjugate

PE. Ex:565nm. Em:574nm.

## Recommended Dilutions

FC 10 µL/Test

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Background

The protein encoded by this gene is the CD3-epsilon polypeptide, which together with CD3-gamma, -delta and -zeta, and the T-cell receptor alpha/beta and gamma/delta heterodimers, forms the T-cell receptor-CD3 complex. This complex plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. The genes encoding the epsilon, gamma and delta polypeptides are located in the same cluster on chromosome 11. The epsilon polypeptide plays an essential role in T-cell development. Defects in this gene cause immunodeficiency. This gene has also been linked to a susceptibility to type I diabetes in women.

## Immunogen Information

### Gene ID

916

### Swiss Prot

P07766

### Immunogen

Recombinant protein of human CD3.

### Synonyms

T3E; TCRE; IMD18; CD3epsilon

## Product Information

### Source

Mouse

### Isotype

IgG1


### Purification

Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.1% Proclin300,0.2% BSA,pH7.2.



Flow cytometry:1X10<sup>6</sup> Human PBMC were surface-stained with PE Mouse anti-Human CD3 mAb (A22795, orange line) or PE Mouse IgG1 isotype control mAb (A22878, blue line).