

A23204

Leader in Biomolecular Solutions for Life Science



RPL22L1 Rabbit pAb

Catalog No.: A23204

Basic Information

Observed MW

Refer to figures

Calculated MW

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable RNA binding activity. Predicted to be a structural constituent of ribosome. Predicted to be involved in cytoplasmic translation. Predicted to be located in ribosome.

Recommended Dilutions

WB 1:1000 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

200916

Swiss Prot

Immunogen

A synthetic peptide corresponding to RPL22L1.

Synonyms

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.