

A23772

Leader in Biomolecular Solutions for Life Science



# EV71 3C Rabbit mAb

Catalog No.: A23772 **Recombinant**

## Basic Information

---

### Observed MW

40kDa

### Calculated MW

20KDa

### Category

SMab Recombinant Monoclonal  
Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Enterovirus

### CloneNo number

ARC59088

## Background

---

## Recommended Dilutions

---

WB 1:10000 - 1:170000

## Immunogen Information

---

Gene ID Swiss Prot

### Immunogen

Recombinant protein of enterovirus EV71 3C.

### Synonyms

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---


Source	Isotype	Purification
Rabbit	IgG	Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.05% proclin300,0.05% BSA,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis lysates of recombinant EV71 3C Protein, using EV71 3C Rabbit mAb (A23772) at 1:160000 dilution.

Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 5ng per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 20s.