

A24476

Leader in Biomolecular Solutions for Life Science



GSDMD (Full Length+N terminal) Rabbit pAb

Catalog No.: A24476

16 Publications

Basic Information

Observed MW

53 kDa/30-32 kDa

Calculated MW

53kDa/32kDa

Category

Polyclonal Antibody

Applications

WB,IF/ICC,IP,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

Gasdermin D is a member of the gasdermin family. Members of this family appear to play a role in regulation of epithelial proliferation. Gasdermin D has been suggested to act as a tumor suppressor. Alternatively spliced transcript variants have been described.

Recommended Dilutions

WB 1:1000 - 1:5000

IF/ICC 1:50 - 1:200

IP 0.5µg-4µg antibody for
200µg-400µg extracts of
whole cells

ELISA Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on
your specific assay
requirements.

Contact



www.abclonal.com

Immunogen Information

Gene ID

69146

Swiss Prot

Q9D8T2

Immunogen

This information is considered to be commercially sensitive.

Synonyms

GSDMD; DF5L; DFNA5L; FKSG10; GSDMDC1; gasdermin-D; GSDMD (Full Length+N terminal)

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of lysates from wild type (WT) and 293T cells transfected with GSDMD (Full Length+N terminal) using GSDMD (Full Length+N terminal) Rabbit pAb (A24476) at 1:1000 dilution incubated overnight at 4°C.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 20 µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020)
.Exposure time: 10s.