

A25864

Leader in Biomolecular Solutions for Life Science



ABflo® 594 Rabbit anti-Mouse CD226/DNAM-1 mAb

Catalog No.: A25864

Basic Information

Observed MW

Calculated MW

38kDa

Category

Monoclonal Antibody

Applications

FC

Cross-Reactivity

Mouse

CloneNo number

ARC67159-ABflo594

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Recommended Dilutions

FC 5 µl per 10⁶ cells in
100 µl volume

Background

Predicted to enable integrin binding activity and protein kinase binding activity. Acts upstream of or within positive regulation of interferon-gamma production and positive regulation of natural killer cell cytokine production. Located in external side of plasma membrane. Is expressed in embryo. Orthologous to human CD226 (CD226 molecule).

Immunogen Information

Gene ID

225825

Swiss Prot

Q8K4F0

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Pta1; DNAM1; DNAM-1; TLISA1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at 2-8°C. Avoid freeze.

Buffer: PBS containing 0.2% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.


Validation Data



Flow cytometry: 1×10^6 293T cells (negative control, left) and 293T (Transfection, right) cells were surface-stained with ABflo® 594 Rabbit anti-Mouse CD226/DNAM-1 mAb (A25864, 5 μ l/Test, orange line) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, blue line). Non-fluorescently stained cells were used as blank control (red line).



Flow cytometry: 1×10^6 293T (Transfection) cells were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Mouse CD226/DNAM-1 mAb (A25864, 5 μ l/Test, right).



Flow cytometry: 1×10^6 C57BL/6 Splenocytes cells were surface-stained with ABflo® 488 Rabbit anti-Mouse CD8a mAb (A23903, 5 μ l/Test) and ABflo® 594 Rabbit anti-Mouse CD226/DNAM-1 mAb (A25864, 5 μ l/Test, left) or ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, right).