

A26096

Leader in Biomolecular Solutions for Life Science



CD5 (Fab) Rabbit mAb

Catalog No.: A26096

Recombinant

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

55kDa

Category

SMab Recombinant Monoclonal
Antibody

Applications

ELISA

Cross-Reactivity

Mouse

Background

This gene encodes a member of the scavenger receptor cysteine-rich (SRCR) superfamily. Members of this family are secreted or membrane-anchored proteins mainly found in cells associated with the immune system. This protein is a type-I transmembrane glycoprotein found on the surface of thymocytes, T lymphocytes and a subset of B lymphocytes. The encoded protein contains three SRCR domains and may act as a receptor to regulate T-cell proliferation. Alternative splicing results in multiple transcript variants encoding different isoforms.

Recommended Dilutions

ELISA

Recommended starting
concentration is 1 µg/mL.
Please optimize the
concentration based on
your specific assay
requirements.

Immunogen Information

Gene ID

12507

Swiss Prot

P06127

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

T1; LEU1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

Fab

Purification

Ni-NTA purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide, 2mM DTT, 0.05% BSA, 50% glycerol, pH7.2.