

A26412-PM

Leader in Biomolecular Solutions for Life Science



ICAM-1/CD54 Rabbit PolymAb®

Catalog No.: A26412-PM

Basic Information

Observed MW

82-100kDa

Calculated MW

58kDa/59kDa

Category

SMab Recombinant Monoclonal
Antibody

Applications

WB,IHC-P,IF/ICC,ELISA

Cross-Reactivity

Human,Mouse,Rat

Background

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

Recommended Dilutions

WB	1:1000 - 1:30000
IHC-P	1:2000 - 1:10000
IF/ICC	1:2000 - 1:4000
ELISA	Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Contact



www.abclonal.com

Immunogen Information

Gene ID

3383/15894

Swiss Prot

P05362/P13597

Immunogen

Recombinant protein of human/mouse ICAM-1/CD54.

Synonyms

MAL; CD54; Icam; Ly-4; Ly-47; Icam-1; MALA-2

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide,0.05% BSA,50% glycerol,pH7.3.