

A28953

Leader in Biomolecular Solutions for Life Science



Anti-Monkey F12 pAb

Catalog No.: A28953

Basic Information

Observed MW

Calculated MW

68 kDa

Category

Polyclonal Antibody

Applications

IF-P

Cross-Reactivity

Cynomolgus monkey

Background

This gene encodes coagulation factor XII which circulates in blood as a zymogen. This single chain zymogen is converted to a two-chain serine protease with an heavy chain (alpha-factor XIIa) and a light chain. The heavy chain contains two fibronectin-type domains, two epidermal growth factor (EGF)-like domains, a kringle domain and a proline-rich domain, whereas the light chain contains only a catalytic domain. On activation, further cleavages takes place in the heavy chain, resulting in the production of beta-factor XIIa light chain and the alpha-factor XIIa light chain becomes beta-factor XIIa heavy chain. Prekallikrein is cleaved by factor XII to form kallikrein, which then cleaves factor XII first to alpha-factor XIIa and then to beta-factor XIIa. The active factor XIIa participates in the initiation of blood coagulation, fibrinolysis, and the generation of bradykinin and angiotensin. It activates coagulation factors VII and XI. Defects in this gene do not cause any clinical symptoms and the sole effect is that whole-blood clotting time is prolonged.

Recommended Dilutions

IF-P 1:200 - 1:400

Immunogen Information

Gene ID

102125559

Swiss Prot

P00748

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

HAF; HAE3; HAEX

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at -20°C

Validation Data



Immunofluorescence analysis of paraffin-embedded Cynomolgus monkey liver tissue using Anti-Monkey F12 pAb (A28953) at a dilution of 1:200 (40x lens).

Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining. High pressure antigen retrieval performed with 0.01M Citrate Buffer (pH 6.0) prior to IF staining.