

A2910

Leader in Biomolecular Solutions for Life Science



EGR4 Rabbit pAb

Catalog No.: A2910

2 Publications

Basic Information

Observed MW

62kDa

Calculated MW

62kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables DNA-binding transcription activator activity, RNA polymerase II-specific and sequence-specific double-stranded DNA binding activity. Involved in positive regulation of transcription by RNA polymerase II. Predicted to be located in nucleoplasm. Predicted to be part of chromatin.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID
1961

Swiss Prot
Q05215

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

AT133; NGFIC; NGFI-C; PAT133; EGR4

Contact



www.abclonal.com

Product Information

Source
Rabbit


Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C

Validation Data



Western blot analysis of lysates from mouse brain, using EGR4 Rabbit pAb (A2910).
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.