

A29244

Leader in Biomolecular Solutions for Life Science



ABflo® 594 Rabbit anti-Mouse CD19 mAb

Catalog No.: A29244

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibody

Applications

FC

Cross-Reactivity

Mouse

Conjugate

ABflo® 594. Ex:588nm. Em:604nm.

Background

Involved in several processes, including B-1 B cell differentiation; positive regulation of phosphatidylinositol 3-kinase activity; and positive regulation of release of sequestered calcium ion into cytosol. Acts upstream of or within B cell receptor signaling pathway. Located in external side of plasma membrane. Is integral component of plasma membrane. Is expressed in liver and spleen. Human ortholog(s) of this gene implicated in common variable immunodeficiency. Orthologous to human CD19 (CD19 molecule).

Recommended Dilutions

FC $\leq 0.25 \mu\text{g}$ per million cells
in 100 μl volume

Immunogen Information

Gene ID
12478

Swiss Prot
P25918

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Contact

 www.abclonal.com

Product Information


Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage
Store at -20°C

Validation Data



Flow cytometry: 1×10^6 C57BL/6 mouse splenocytes were surface-stained with ABflo® 594 Rabbit IgG isotype control (A23821, 5 μ l/Test, left) or ABflo® 594 Rabbit anti-Mouse CD19 mAb (A29244, 0.25 μ g, right).