

A29344

Leader in Biomolecular Solutions for Life Science



PE Rabbit anti-Mouse CD8a mAb

Catalog No.: A29344

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibody

Applications

FC

Cross-Reactivity

Mouse

Conjugate

PE. Ex:565nm. Em:574nm.

Background

Enables identical protein binding activity. Acts upstream of or within several processes, including cytotoxic T cell differentiation; defense response to virus; and positive regulation of calcium-mediated signaling. Located in external side of plasma membrane. Is expressed in several structures, including endocrine gland; exocrine gland; genitourinary system; gut; and hemolymphoid system gland. Orthologous to human CD8A (CD8a molecule).

Recommended Dilutions

FC $\leq 0.25 \mu\text{g}$ per million cells
 in 100 μl volume

Immunogen Information

Gene ID
12525

Swiss Prot
P01731

Immunogen

Recombinant protein (or fragment). This information is considered to be commercially sensitive.

Synonyms

Ly-2; Ly-B; Ly-35; Lyt-2

Contact

 www.abclonal.com

Product Information

Source
Rabbit


Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C

Validation Data



Flow cytometry: 1X10⁶ C57BL/6 mouse splenocytes were surface-stained with PE Rabbit IgG isotype control (A24172,5 µl/Test,left) or PE Rabbit anti-Mouse CD8a mAb (A29344,0.25 µg,right).