

A2942

Leader in Biomolecular Solutions for Life Science



# CMTM5 Rabbit pAb

Catalog No.: A2942

## Basic Information

---

### Observed MW

Refer to figures

### Calculated MW

25kDa

### Category

Polyclonal Antibody

### Applications

### Cross-Reactivity

Human

## Background

---

This gene encodes a member of the chemokine-like factor superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. The encoded protein may exhibit tumor suppressor activity. Alternative splicing results in multiple transcript variants.

## Recommended Dilutions

---

## Immunogen Information

---

### Gene ID

116173

### Swiss Prot

Q96DZ9

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

CKLFSF5; CMTM5

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.