

A3001

Leader in Biomolecular Solutions for Life Science



BPIFB1 Rabbit pAb

Catalog No.: A3001

Basic Information

Observed MW

Refer to figures

Calculated MW

52kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

The protein encoded by this gene may be involved in the innate immune response to bacterial exposure in the mouth, nasal cavities, and lungs. The encoded protein is secreted and is a member of the BPI/LBP/PLUNC protein superfamily. This gene is found with other members of the superfamily in a cluster on chromosome 20.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

92747

Swiss Prot

Q8TDL5

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

LPLUNC1; C20orf114; BPIFB1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.