

A3574

Leader in Biomolecular Solutions for Life Science



# CCDC158 Rabbit pAb

Catalog No.: A3574

## Basic Information

### Observed MW

135kDa

### Calculated MW

127kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Mouse

## Background

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

339965

### Swiss Prot

Q5M9N0

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

CCDC158

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of lysates from mouse testis, using CCDC158 Rabbit pAb (A3574).  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.