

A3665

Leader in Biomolecular Solutions for Life Science



MORF4L1 Rabbit pAb

Catalog No.: A3665

Basic Information

Observed MW

41kDa

Calculated MW

41kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables protein N-terminus binding activity. Involved in double-strand break repair via homologous recombination and histone modification. Located in nuclear speck. Part of NuA4 histone acetyltransferase complex and Sin3 complex.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

10933

Swiss Prot

Q9UBU8

Immunogen

A synthetic peptide of human MORF4L1

Synonyms

Eaf3; MEAF3; MRG15; FWP006; S863-6; HsT17725; MORFRG15; MORF4L1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.