

A4330

Leader in Biomolecular Solutions for Life Science



SAFB2 Rabbit pAb

Catalog No.: A4330

Basic Information

Observed MW

160kDa

Calculated MW

107kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene, along with its paralog (scaffold attachment factor B1), is a repressor of estrogen receptor alpha. The encoded protein binds scaffold/matrix attachment region (S/MAR) DNA and is involved in cell cycle regulation, apoptosis, differentiation, the stress response, and regulation of immune genes.

Recommended Dilutions

WB 1:1000 - 1:2000

Immunogen Information

Gene ID

9667

Swiss Prot

Q14151

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-255 of human SAFB2 (NP_055464.1).

Synonyms

SAFB2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data

