

A4605

Leader in Biomolecular Solutions for Life Science



RAB11FIP5 Rabbit pAb

Catalog No.: A4605

Basic Information

Observed MW

65kDa

Calculated MW

70kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

Enables gamma-tubulin binding activity. Involved in cellular response to acidic pH; negative regulation of adiponectin secretion; and regulation of protein localization to cell surface. Located in centriolar satellite and mitochondrial outer membrane.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

26056

Swiss Prot

Q9BXF6

Immunogen

A synthetic Peptide of human RAB11FIP5

Synonyms

GAF1; pp75; RIP11; gaf-1; rab11-FIP5; RAB11FIP5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG


Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of extracts of HeLa cells, using RAB11FIP5 antibody (A4605).
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.