

A4786

Leader in Biomolecular Solutions for Life Science



# RHOF Rabbit pAb

Catalog No.: A4786

## Basic Information

---

### Observed MW

27kDa

### Calculated MW

24kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse

## Background

---

Predicted to enable GTP binding activity; GTPase activity; and protein kinase binding activity. Involved in actin filament organization. Located in extracellular exosome.

## Recommended Dilutions

---

WB 1:500 - 1:2000

## Immunogen Information

---

### Gene ID

54509

### Swiss Prot

Q9HBH0

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 40-120 of human RHOF (NP\_061907.2).

### Synonyms

RIF; ARHF; RHOF

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of various lysates using RHO Rabbit pAb (A4786) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 90s.