

A4901

Leader in Biomolecular Solutions for Life Science



# ZNF410 Rabbit pAb

Catalog No.: A4901

## Basic Information

### Observed MW

57kDa

### Calculated MW

52kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

Enables sequence-specific double-stranded DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

## Recommended Dilutions

WB 1:500 - 1:1000

## Immunogen Information

### Gene ID

57862

### Swiss Prot

Q86VK4

### Immunogen

A synthetic Peptide of human ZNF410

### Synonyms

APA1; APA-1; ZNF410

## Contact

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification


Affinity purification

### Storage

Store at 4°C. Avoid freeze / thaw cycles.  
Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of extracts of HeLa cells, using ZNF410 antibody (A4901).  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.