

A4904

Leader in Biomolecular Solutions for Life Science



# PCTP Rabbit pAb

Catalog No.: A4904

## Basic Information

---

### Observed MW

25kDa

### Calculated MW

25kDa

### Category

Polyclonal Antibody

### Applications

WB

### Cross-Reactivity

Human

## Background

---

Enables phosphatidylcholine binding activity and phosphatidylcholine transporter activity. Involved in phospholipid transport. Predicted to be located in cytosol.

## Recommended Dilutions

---

WB 1:500 - 1:1000

## Immunogen Information

---

### Gene ID

58488

### Swiss Prot

Q9UKL6

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

PC-TP; STARD2; PCTP

## Contact

---

 [www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

## Validation Data

---



Western blot analysis of lysates from HL-60 cells, using PCTP Rabbit pAb (A4904).  
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.