

A4951

Leader in Biomolecular Solutions for Life Science



UBE2Z Rabbit pAb

Catalog No.: A4951

Basic Information

Observed MW

38kDa

Calculated MW

38kDa

Category

Polyclonal Antibody

Applications

WB

Cross-Reactivity

Human

Background

This gene encodes an enzyme which ubiquitinates proteins which participate in signaling pathways and apoptosis.

Recommended Dilutions

WB 1:500 - 1:1000

Immunogen Information

Gene ID

65264

Swiss Prot

Q9H832

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

USE1; HOYS7; UBE2Z

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, pH7.3.

Validation Data



Western blot analysis of lysates from MDA-MB-453 cells, using UBE2Z Rabbit pAb (A4951).

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 μ g per lane.

Blocking buffer: 3% nonfat dry milk in TBST.