

A4969

Leader in Biomolecular Solutions for Life Science



# ADIPOR2 Rabbit pAb

Catalog No.: A4969

## Basic Information

### Observed MW

### Calculated MW

44kDa

### Category

Polyclonal Antibody

### Applications

### Cross-Reactivity

Human, Mouse

## Background

The adiponectin receptors, ADIPOR1 (MIM 607945) and ADIPOR2, serve as receptors for globular and full-length adiponectin (MIM 605441) and mediate increased AMPK (see MIM 602739) and PPAR-alpha (PPARA; MIM 170998) ligand activities, as well as fatty acid oxidation and glucose uptake by adiponectin (Yamauchi et al., 2003 [PubMed 12802337]).

## Recommended Dilutions

## Immunogen Information

### Gene ID

79602

### Swiss Prot

Q86V24

### Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

### Synonyms

PAQR2; ACDCR2; ADIPOR2

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.